

CareVis

Single-use elastomeric infusion pumps for
pain therapy



PAJUNK®


PROMECON

CareVis elastomeric pumps

Models

Infusion pumps for a variety of indications

Elastomeric infusion pumps for single use

CareVis® is an elastomeric pump for single use for the infusion of cytostatics, antibiotics, pain medication and more. The pump offers patients very high mobility during treatment and is a safe and simple alternative to electronic devices.

CareVis® can be used in both stationary and ambulatory settings and is available with continuous flow rates between 5 ml/hr and 10 ml/hr and in sizes between 275 ml and 600 ml. In addition, there are PCA models with bolus button for the patient and bolus sizes up to 2,0 ml as well as different lock-out times and multirate models with removable adjustment key.



Advantages

- easy application, works without power source
- large range of infusion rates, filling volumes and special models, optimally adapted to treatment needs
- safety and easy handling due to closed single use system
- shockproof and hard case to protect the balloon and exclude external influences
- high delivery accuracy
- air and particle filter as well as king resistant cable ensure maximum safety
- bacteria proof vent cap - priming with closed cap without dripping
- enable ambulatory treatment and early discharge
- cost-effective
- transparent balloon and case for clear visual control of the medication solution
- minimal dead space for complete drug administration
- easy to read volume indicator
- free of latex and DEHP

Indication

Short, medium, or longterm infusions for pain medications.

CareVis® is suitable for intravenous, subcutaneous, epidural and tissue infusion.

Suitable for regional anaesthesia, epidural, on nerves, and intralesionally in surgical wounds. The CareVis® is also available with neuraxial connector according to ISO 80369-6.

Models


Infusion pumps for a variety of indications

CareVis - continuous basal flow

Infusion pump with predefined flow rate, universally applicable for continuous applications. Light weight and easy to use.

Application:

- Wide range of applications for continuous infusion including keeping access open, intravenous, subcutaneous, epidural, into tissue for nerve blocks and wound infiltration.



CareVis PCA - continuous basal rate and bolus function

With continuous basal rate and bolus button for patient controlled analgesia.

Application:

- Postoperative pain therapy, systemic or regional



CareVis BoluS - bolus funktion only

Medication is supplied exclusively via the bolus function without a continuous basal rate

Application:

- Postoperative pain therapy, systemic or regional



CareVis VariO - seven adjustable flow rates

Multirate version with seven adjustable flow rates for basal flow selection.

Application:

- full spectrum of regional anaesthesia including wound infiltration.



CareVis VariO Plus - adjustable flow rate inclusive bolus function

Combines the seven adjustable flow rates of the CareVis VariO multirate function with a bolus function.

Application:

- Regional anaesthesia



Also available in **NRFit**

CareVis elastomeric pumps

Pumps

CareVis pumps

Description	Item No.	NRFit Item No.	PU
450ml (315ml /480ml) basal rate 5ml/h preset, cont. basal rate	X50	X50-N	10
450ml (315ml /480ml) basal rate 8ml/h preset, cont. basal rate	X80	X80-N	10
450ml (315ml /480ml) basal rate 10ml/h preset, cont. basal rate	X100	X100-N	10
600ml (420ml /650ml) basal rate 8ml/h preset, cont. basal rate	G80	G80-N	10
600ml (420ml /650ml) basal rate 10ml/h preset, cont. basal rate	G100	G100-N	10

CareVis PCA pumps

Description	Item No.	NRFit Item No.	PU
450ml (315ml /480ml) basal rate 5ml/h bolus: 2 ml, LO 15 min preset, cont. basal rate	X5015W	X5015W-N	10
450ml (315ml /480ml) basal rate 8ml/h bolus: 2 ml, LO 15 min preset, cont. basal rate	X8015W	X8015W-N	10

CareVis BoluS pumps

Description	Item No.	NRFit Item No.	PU
275ml (190ml /310ml) bolus: 2 ml, LO 30 min	L0030W	L0030W-N	10

CareVis VariO pumps

Description	Item No.	NRFit Item No.	PU
450ml (315ml /480ml) basal rate 0-2-4-6-8-10-12-14 ml/h adjustable, cont. basal rate	VX214	VX214-N	10



Also available in NRFit®

CareVis VariO Plus pumps

Description	Item No.	NRFit Item No.	PU
450ml (315ml /480ml) basal rate 0-2-4-6-8-10-12-14 ml/h bolus: 2 ml, LO 15 min adjustable, cont. basal rate	VX21415W	VX21415W-N	10
450ml (315ml /480ml) basal rate 0-2-4-6-8-10-12-14 ml/h bolus: 2 ml, LO 30 min adjustable, cont. basal rate	VX21430W	VX21430W-N	10

Pump bags

Description	Item No.	PU
Carrying bag for 450-600 ml pumps	TT003	30



Also available in NRFit®

NOTES

A series of horizontal dashed lines for taking notes.

NOTES

A series of 25 horizontal dashed lines for writing notes.



Listed elastomeric pumps of Promecon exclusively available from PAJUNK®

Please contact us for any questions at phone +49 7704 80008-0 or by e-mail at info@pajunk-medipro.com.

Did we catch your interest?
You can find more information about Pajunk products at: www.pajunk.eu

Promecon is medical device manufacturer with a certified quality management system (ISO 13485:2016).

This catalogue is intended solely for information about our products. Errors excepted. Catalogue no. XS400144EU

PAJUNK®

Distribution company for:
Germany, Austria, BeNeLux & Poland

PAJUNK® Medical Produkte GmbH
Central Europe

Karl-Hall-Straße 1
78187 Geisingen, Germany
Tel.: +49 7704 8008-0
Fax: +49 7704 8008-150
info@pajunk-medipro.com
pajunk.de & pajunk.eu

Distribution company for:
Switzerland

PAJUNK® Swiss AG
Schwitzerland

Neuhofstrasse 5A
6340 Baar, Schwitzerland
Tel.: +41 41 760 96 27

info@pajunk.ch
pajunk.ch



PROMECON

PROMECON GmbH | Kreuzerstraße 10 | 20095 Hamburg
Tel.: +49 40 3690 1690 | Fax: +49 40 3690 1699
info@promecon-medical.de | www.promecon-medical.com

NRFit® is a registered trademark of GEDSA used with their permission.

CareVis® is a registered trademark of Promecon used with their permission.

Certification Experts B.V.
Amerlandsdweg 7, 3621 ZC Breukelen,
The Netherlands



CE 1370

S&S Med Co., Ltd.
67, Bongseong-ro, Gunpo-si
Gyeonggi-do, 15848, Korea



XS400144EU 05/22