

**PAJUNK®**

## ***MultiTip***

*Economic instrument system  
with Multi inserts*



Laparoscopic instruments

MADE IN GERMANY

## Instrument inserts MultiTip

# The perfect combination of cost-effectiveness and sustainability

PAJUNK® is the right choice for anyone who wishes to avoid the amount of waste from disposable inserts but who still wants a cost-effective alternative to reusable inserts. With MultiTip PAJUNK® offers a cost-effective solution, which can be reprocessed for five additional uses. The corresponding jaw inserts are delivered sterile for immediate use.


- ➔ Secure connection system
- ➔ Instrument change in the OR
- ➔ No regrinding
- ➔ Cost-effective and environmentally friendly



Multi scissors inserts, MultiTip




Multi grasping forceps inserts, MultiTip




Guide tube

Ø 5 mm



Drawbar

Ø 5 mm



A combined click and screw system helps avoid the undesired loosening of the inserts.

### Instrument system with Multi inserts



#### Multi scissors inserts, Ø 5 mm

Model

Item.-No.

PU

MultiTip METZENBAUM L

double-action, jaw insert length: 20 mm

large

1293-94-01

10



MultiTip METZENBAUM M

double-action, jaw insert length: 15,5 mm

medium

1293-94-02

10



MultiTip METZENBAUM S

double-action, jaw insert length: 10 mm

small

1293-94-03

10



MultiTip METZENBAUM P

double-action, jaw insert length: 16,5 mm

medium/  
pointed

1293-94-04

10



MultiTip DISSECT M

double-action, jaw insert length: 15,5 mm

medium

1293-94-12

10



#### Multi grasping forceps inserts, Ø 5 mm

Item.-No.

PU

MultiTip ENDOCLINCH

double-action, jaw insert length: 21 mm

1293-54-10

10



MultiTip MARYLAND

double-action, jaw insert length: 20 mm

1293-54-20

10



ErgoSys handle  
ErgoFlex

#### ErgoSys handle

Item.-No.

PU

ErgoFlex handle for modular instrument inserts

1292-30-00

1

#### Guide tubes, Ø 5 mm

Working length

Item.-No.

PU

TipTube

330 mm

1293-15-05

1

TipTube

450 mm

1293-16-05

1

#### Drawbars, Ø 5 mm

Working length

Item.-No.

PU

TipRod Ergo

330mm

1293-20-05

1

TipRod Ergo

450mm

1293-21-05

1

#### PAJUNK® GmbH

Medizintechnologie

Karl-Hall-Strasse 1

D-78187 Geisingen/Germany

Phone +49 (0) 77 04/92 91-0

Telefax +49 (0) 77 04/92 91-6 00

www.pajunk.com

#### PAJUNK® Medical Produkte GmbH

D.A.CH • BeNeLux

Karl-Hall-Strasse 1

D-78187 Geisingen/Germany

Phone +49 (0) 77 04/80 08-0

Telefax +49 (0) 77 04/80 08-150

www.pajunk.com