


PAJUNK®



REGIONAL ANAESTHESIA
& PAIN MANAGEMENT

UniPlex®

*Stimulating Single Shot
Nerve Block Needles*



The Stimulating Needle with NanoLine® Coating

The UniPlex NanoLine needles have been developed for single shot nerve block applications, allowing for stimulation guidance only. With the introduction of the NanoLine thin layer coating technology, PAJUNK® became the innovative technology leader in stimulation guided nerve block needles. The new type of coating offers improved gliding properties and precise stimulation while reducing needle thickness.

→ All coated Pajunk nerve block needles use NanoLine technology offering superior gliding properties.

i

AVAILABLE OPTIONS

UniPlex II NRFit®

Stimulating single shot nerve block needles with non-removable injection tubing, removable stimulation cable, colour-coded needle hub and NanoLine Coating

More information at: pajunk.com

NanoLine Coating

Very thin polymer film, smooth surface, accurate inner and outer diameter

- Excellent puncture and gliding properties through smooth surface
- Improves visibility under ultrasound¹
- Precise stimulation through the non-insulated needle tip

Depth Markings

- Easy to read and identify

Echogenic Needle Tip

Facet tip with two inclination angles

- Improves needle tip visibility under ultrasound



Also available in NRFit®

UniPlex®

UniPlex needle with facet tip

Size	UniPlex Item no.	UniPlex II NRFit® Item no.	PU
25G x 35 mm (1 3/8")	001156-80	001246-80	10
25G x 50 mm (2")	001156-81	001246-81	10
25G x 80 mm (3 1/5")	001156-86		10
24G x 25 mm (1")	001156-75	001246-75*	10
24G x 50 mm (2")	001156-85		10
22G x 40 mm (1 5/8")	001156-70	001246-70	10
22G x 50 mm (2")	001156-74	001246-74	10
22G x 80 mm (3 1/5")	001156-71	001246-71	10
22G x 100 mm (4")	001156-84		10
22G x 120 mm (4 3/4")	001156-82	001246-82	10
21G x 100 mm (4")	001156-77	001246-77	10
20G x 120 mm (4 3/4")	001156-72	001246-72	10
20G x 150 mm (6")	001156-76	001246-76	10
20G x 180 mm (7 1/4")	001156-73		10

UniPlex needle with facet S tip

Size	UniPlex Item no.	UniPlex II NRFit® Item no.	PU
25G x 35 mm (1 3/8")	001157-80		10
22G x 40 mm (1 5/8")	001157-70		10
22G x 50 mm (2")	001157-74	001247-74	10
21G x 100 mm (4")	001157-77	001247-77	10

UniPlex needle with SPROTTE® tip

Size	UniPlex Item no.	UniPlex II NRFit® Item no.	PU
25G x 50 mm (2")	001156-29G		10
24G x 40 mm (1 5/8")	001156-30G		10
22G x 50 mm (2")	001156-31G	001246-31G	10
22G x 70 mm (2 3/4")	001156-31H		10
22G x 90 mm (3 1/2")	001156-31J	001246-31J	10
22G x 150 mm (6")	001156-28L	001246-28L	10

*24G x 30 mm (1 1/5")