

**PAJUNK®**

## ***Customized Procedural Kits***


*PAJUNK® customized kits -  
create your own kit*



Customized kits

## Customized kits

### ***In 5 steps to your own customized kit***



#### ***Why a customized kit?***

Not only patient safety but also time plays an increasingly important role in everyday clinical practice. In order to optimize and perfect these processes, all components of a process chain have to be optimally coordinated with each other. With the PAJUNK® customized procedural kits, you do not only combine high quality in products, but also optimal process improvement in the respective areas:

- Spinal anesthesia
- Epidural anesthesia
- Peripheral nerve blocks
- Pain therapy

The configuration of a PAJUNK® customized kit is performed in a modular system. According to the clinical procedure, a basic procedural framework is assembled from a large number of building blocks. Building block for building block is then equipped with the appropriate products. This kind of modular system creates a variety of combination possibilities for your customized kit. With this system, PAJUNK® adapts to your requirements and needs in the medical field. Based on the extensive PAJUNK® portfolio, we configure your kit.

## ***Customized kits at PAJUNK® - Create your own kit***

## Customized kits

# PAJUNK® - Service all along the line

In order to meet your desires and requirements for a customized kit, our staff members are at your disposal.



### Services at a glance:

- Competent service and advice from our staff members
- Presentation of the relevant products at the first appointment
- Free, sterile or non-sterile samples for clinical testing or inspection.
- Consideration of individual customer or component requirements

*The service makes  
the difference*

- ➡ Individual solutions – the right kit for everyone
- ➡ All components for the procedure at your disposal
- ➡ Safety due to standardization

## Customized kits

# NRFit connectors for more patient safety

NRFit™




Due to various incidents in the past, which were based on misconnection of different patient lines, the ISO standard 80369-6 has been adopted in 2016. It defines new requirements for connectors in the area of regional anesthesia.

This standard is intended to reduce the risk of accidental misconnections and injection, thus increasing patient safety.

For this reason, PAJUNK® also offers the complete portfolio with NRFit connectors.


Overview of portfolio expansion on NRFit connectors:



- ➔ Regional anesthesia components become NRFit
- ➔ Local anesthesia components remain Luer

**To increase the process reliability – especially in the early conversion phase – we recommend the assembly of a customized kit in which both connector types are combined.**

## Accessories for your customized kit




Preparation	Puncture	Postprocessing
<ul style="list-style-type: none"><li>■ Protective clothing</li><li>■ Disinfection and hygienics</li><li>■ Patient cover</li><li>■ Injection accessories (syringes &amp; cannulas)</li><li>■ Filter cannulas &amp; filter straws</li><li>■ Ultrasound probe covers &amp; ultrasound gel</li><li>■ Swab cups &amp; disinfection cups</li></ul>	<ul style="list-style-type: none"><li>■ <b>Spinal anesthesia:</b> SPROTTE® cannula IntraLong</li><li>■ <b>Epidural anesthesia:</b> EpiLong EpiSpin Lock</li><li>■ <b>Peripheral nerve block:</b> SonoSystem Sono cannulas E-Cath acc. Tsui StimuLong Sono Plexus NanoLine</li><li>■ <b>Pain therapy:</b> InfiltraLong &amp; FuserPump Thermolesion cannulas</li></ul>	<ul style="list-style-type: none"><li>■ Catheter fixation</li><li>■ Surgical sewing material</li><li>■ Wound dressing &amp; plasters</li></ul>

## Preparation

In order to reduce the risk of contamination, increase the patient's safety and meet the requirements for sterility some preparations must be taken. The following building blocks are located in the PAJUNK® portfolio:

### Disinfection

- 
- ➔ Round swabs
  - ➔ Cotton swab tweezers and forceps
  - ➔ Swabsticks
  - ➔ Compresses


### Drapes

- ➔ Blue, green or transparent fenestrated drapes
- ➔ Blue, green or transparent drapes
- ➔ Fenestrated drapes with transparent window and adhesive stripe
- ➔ Fenestrated drapes with special forms

### Sterile work

- ➔ Medical gloves
- ➔ Surgical gowns


### Injection

- 
- ➔ Injection cannulas & safety cannulas
  - ➔ Syringes Luer Lock /Slip & NRFit Lock /Slip
  - ➔ Filter cannulas & filter straws Luer / NRFit


## Puncture – Kit components for each indication

The optimal choice of puncture equipment is indispensable for every clinician. Not only due to individual preferences, but also in terms of quality and quantity. PAJUNK® offers the following portfolio for your customized kit.

### Spinal anesthesia

- 
- ➔ SPROTTE® spinal cannula
  - ➔ IntraLong

### Epidural anesthesia

- 
- ➔ EpiLong I
  - ➔ EpiLong Tuohy
  - ➔ EpiLong Soft
  - ➔ EpiLong Soft Sono
  - ➔ Tuohy cannula
  - ➔ EpiSpin Lock

## *Punction – Kit components for each indication*

### Peripheral nerve blocks

- ➔ PlexoLong NanoLine
- ➔ SonoSystem
- ➔ Sono cannulas
- ➔ E-Cath acc. Tsui



### Pain therapy

- ➔ RF thermolesion cannulas
- ➔ RF Sono thermolesion cannulas
- ➔ X-ray contrast cannulas
- ➔ InfiltraLong
- ➔ FuserPump



### *Post processing*

A wide range of components is available for catheter fixation and wound care. From conventional plaster over wound dressings to suture material you can find everything in our portfolio.



### Catheter fixation & wound care

- ➔ FixoCath catheter fixation
- ➔ Plaster stripes
- ➔ Transparent film dressings
- ➔ Suture material





## Customized kits at PAJUNK® - Create your own kit



Your own  
customized kit

**PAJUNK® GmbH**  
Medizintechnologie  
Karl-Hall-Straße 1  
D-78187 Geisingen  
phone +49 (0) 77 04 /92 91-0  
fax +49 (0) 77 04 /92 91-6 00  
[www.pajunk.com](http://www.pajunk.com)

**PAJUNK® Medical Produkte GmbH**  
D•A•CH • BeNeLux  
Karl-Hall-Straße 1  
D-78187 Geisingen  
phone +49 (0) 77 04 / 80 08-0  
fax +49 (0) 77 04 /80 08-150  
[www.pajunk.com](http://www.pajunk.com)

XS400108 A 09/18